

**IEC Tipo K 9kg Pull Power (herradura) Termopar  
magnético, cable aislado PFA con revestimiento  
trenzado de acero inoxidable que termina en  
miniatura con enchufe estándar**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

## Specifications

<b>General Description</b>	Rapid, reliable temperature sensing on ferrous surfaces
<b>Sensor Type</b>	Type K thermocouple to IEC 584
<b>Thermocouple Type</b>	K
<b>Thermocouple Junction</b>	Grounded at tip
<b>Magnet</b>	25 x 25 x 40mm (HxWxL) with powerful 9kg pull
<b>Cable Termination</b>	1.5 or 3 metres PFA Teflon® insulated with stainless steel over-braid, fitted plug termination
<b>Termination Type</b>	Fitted Mini or Standard Plug
<b>Max. Temperature</b>	+ 250°C
<b>Min. Temperature</b>	-50°C