


Sonda Pt100 de uso general con cable PFA - Tipo RTD / PRT



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Sonda Pt100 de uso general con cable PFA - Tipo RTD / PRT

Sonda multipropósito Pt100 según IEC751

Clase B, conexión de 4 hilos

Funda rígida de acero inoxidable 316, 3,0, 4,0, 4,5 o 6,0 mm de diámetro

1 o 2 metros de cable apantallado y aislado en teflón

Clasificación de temperatura del sensor: de -50 C a +250 C

Características técnicas**Specifications**

General Description	Multipurpose Pt100 Probe to IEC751
Sensor Type	Pt100 (100 Ohms @ 0°C) to IEC 751, Class B, 4 wire
Sheath Material	3.0, 4.0, 4.5 or 6.0mm diameter 316 stainless steel sheath.
Cable Length	1 or 2 metre flexible Teflon® insulated, screened
Termination Type	4-wire tails allows connection to any PT100 instrument 2, 3 or 4 wire
Max. Temperature	+250°C
Min. Temperature	-50°C
Standards Met	IEC751