

## Accesorios de compresin de acero inoxidable - rosca mtrica (serie M)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Para adaptarse a una amplia gama de sensores con aislamiento mineral y otros tipos de sonda ofrecidos por Labfacility.

Los accesorios de compresin de acero inoxidable con roscas mtricas (serie M) se utilizan ampliamente para conectar tuberías y tubos en diversas aplicaciones industriales.

Estos accesorios cuentan con roscas mtricas, estandarizadas para garantizar la compatibilidad con los sistemas de tuberías mtricas.

Fabricados en acero inoxidable, ofrecen una resistencia a la corrosin, durabilidad y resistencia superiores, lo que los hace adecuados para entornos hostiles como el procesamiento químico, los productos farmacuticos y las industrias de alimentos y bebidas.

### Características técnicas

## Specifications

<b>General Description</b>	Thermocouple installation fittings with Metric Thread
<b>Material</b>	Male in-line stainless steel compression type with nut & olive
<b>Thread</b>	M8, M16 & M20 thread available.
<b>Probe Diameter</b>	Suits 1.0mm, 1.5mm, 2.0mm, 3.0mm, 4.5mm, 6.0mm or 8.0mm probe diameters