

Accesorios de compresin de acero inoxidable - rosca mtrica (serie M)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Para adaptarse a una amplia gama de sensores con aislamiento mineral y otros tipos de sonda ofrecidos por Labfacility.

Los accesorios de compresin de acero inoxidable con roscas mtricas (serie M) se utilizan ampliamente para conectar tuberías y tubos en diversas aplicaciones industriales.

Estos accesorios cuentan con roscas mtricas, estandarizadas para garantizar la compatibilidad con los sistemas de tuberías mtricas.

Fabricados en acero inoxidable, ofrecen una resistencia a la corrosin, durabilidad y resistencia superiores, lo que los hace adecuados para entornos hostiles como el procesamiento químico, los productos farmacuticos y las industrias de alimentos y bebidas.

Características técnicas

Specifications

General Description	Thermocouple installation fittings with Metric Thread
Material	Male in-line stainless steel compression type with nut & olive
Thread	M8, M16 & M20 thread available.
Probe Diameter	Suits 1.0mm, 1.5mm, 2.0mm, 3.0mm, 4.5mm, 6.0mm or 8.0mm probe diameters