

Termopar de horno ceramico para alta temperatura (1400 C) con funda de porcelana aluminosa - Tipo N - 300 mm de longitud



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Termopares de horno ceramico para alta temperatura. El tubo de soporte de acero inoxidable de 100 mm de largo fija la funda cermica al cabezal industrial y est incluido en la longitud total de la funda.

Caractersticas tcnicas

Specifications

Product Code	XE-6953-001
General Description	Ceramic thermocouples for furnaces, kilns and other high temperature applications
Sensor Type	Made with rugged 20 gauge (0.8mm dia) thermocouple wire for long-term performance
Thermocouple Type	Type N with Aluminous Porcelain sheaths, rated to 1400°C
Sheath Material	12 mm diameter ceramic sheath
Max. Temperature	+1400°C
IP Rating	IP67 rated industrial protection head with cable entry and screw terminal block
Length	300mm
Support Tube	IP67 rated industrial protection head with cable entry and screw terminal block