

## Sonda de aire mÃfÃ³vil de mano - Tipo K,T



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Cubierta de acero inoxidable con asa moldeada y 2 metros de cable de extensión enrollado terminado en un tapón de termopar en miniatura. Las sondas de aire termopar son ideales para medir la temperatura del aire quieto o en movimiento y los gases.

## Specifications

<b>General Description</b>	Handheld Moving Air Probe, Ideal for the measurement of air and gases
<b>Sensor Type</b>	Funnelled shield allows for maximum airflow across measuring junction
<b>Thermocouple Type</b>	K or T
<b>Calibration</b>	IEC
<b>Probe Diameter</b>	30mm
<b>Probe Length</b>	120mm
<b>Cable Length</b>	2m
<b>Cable Termination</b>	Fitted Plug