

## Sonda de aire mvil de mano - Tipo K IEC 30 mm x 120 mm



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Funda de acero inoxidable con mango moldeado y 2 metros de cable de extensión en espiral terminado en un enchufe de termopar en miniatura. Las sondas de aire de termopar son ideales para medir la temperatura de aire y gases quietos o en movimiento.

### Características técnicas

## Specifications

<b>Product Code</b>	XE-3845-001
<b>General Description</b>	Handheld Moving Air Probe, Ideal for the measurement of air and gases
<b>Sensor Type</b>	Funnelled shield allows for maximum airflow across measuring junction
<b>Thermocouple Type</b>	K IEC
<b>Probe Diameter</b>	30mm
<b>Probe Length</b>	120mm
<b>Cable Length</b>	2m
<b>Cable Termination</b>	Fitted plug