

Sonda de aire mvil de mano tipo T IEC 30 mm de dimetro x 120 mm de longitud



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Funda de acero inoxidable con mango moldeado y 2 metros de cable de extensión en espiral terminado en un enchufe de termopar en miniatura. Las sondas de aire de termopar son ideales para medir la temperatura de aire y gases quietos o en movimiento.

Características técnicas

Specifications

Product Code	XE-3846-001
General Description	Handheld Moving Air Probe, Ideal for the measurement of air and gases. Funnelled shield allows for maximum airflow across measuring junction
Thermocouple Type	T IEC
Probe Diameter	30mm
Probe Length	120mm
Cable Length	2m
Cable Termination	Fitted plug