

## Termopares con aislamiento mineral con cabezal de terminal compacto KNS - Tipo K



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

### **?Necesita su sensor de temperatura calibrado?**

Trazable según los estándares de calibración ISO17025 (UKAS), Labfacility ofrece una opción de calibración trazable de 3 o 5 puntos utilizando nuestra instalación de calibración interna, puede seleccionar cualquier punto de temperatura entre -10 oC y +200 oC.

**Tenga en cuenta que si agrega una opción de calibración a su pedido, esto agregará 2-3 días de tiempo de entrega a su pedido.**

## **Specifications**

<b>General Description</b>	Flexible mineral insulated probe with 310 stainless steel sheath
<b>Sensor Type</b>	Type 'K' (Nickel Chromium/Nickel Aluminium) to IEC 584
<b>Thermocouple Junction</b>	Insulated hot junction
<b>Element/hot junction</b>	Single element, junction insulated from sheath (offers protection against spurious electrical signals)
<b>Probe Diameter</b>	6.0mm diameter
<b>Sheath Material</b>	310 stainless steel sheath
<b>Head Type</b>	KNS, Cast aluminium finish, screw type lid with chain. Head fitted 2-way ceramic block
<b>Gland</b>	M16 x 1.5mm Plated brass gland fitted
<b>Max. Temperature</b>	+1100°C (Max temp of head is 170°C)
<b>Min. Temperature</b>	-40°C
<b>IP Rating</b>	IP67