



Sonda de temperatura fija de los neumticos de competicin



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Specifications

Product Code	XE-3891-001
General Description	Can be used with any digital thermometer with a type k mini input socket
Sensor Type	Fast response fixed racing tyre probe with a max depth of six millimetres
Thermocouple Type	K
Calibration	IEC
Probe Length	max depth of 6mm
Cable Length	1m
Cable Termination	Miniature Plug
Max. Temperature	+250°C
Min. Temperature	-50°C
Features #1	Junction Type: Grounded