

## Carcasa remota montada en pared de 4-20mA, entrada Pt100 con cable PFA o SSOB y sensor PRT



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Carcasa moldeada en ABS de uso general, los puntos de fijación externos proporcionan un medio simple y rápido de montaje de sensores directamente en los conductos, dispensando con la necesidad de placas de montaje separadas.

**?Necesita su sensor de temperatura calibrado?**

Trazable según los estándares de calibración ISO17025 (UKAS), Labfacility ofrece una opción de calibración trazable de 3 o 5 puntos utilizando nuestra instalación de calibración interna, puede seleccionar cualquier punto de temperatura entre -10 oC y +200 oC.

**Tenga en cuenta que si agrega una opción de calibración a su pedido, esto agregará 2-3 días de tiempo de entrega a su pedido.**

## Specifications

<b>General Description</b>	Compatible with Pt100 platinum resistance thermometers
<b>Sensor Type</b>	Surface Mounting for remote installation
<b>Cable Length</b>	1 or 2 metre
<b>Cable Type</b>	PFA screened or Fiberglass stainless steel overbraid available
<b>Core / Strands</b>	7/0.2mm
<b>Head Type</b>	High Impact ABS Plastic Construction
<b>Thread</b>	Entry Threads M16
<b>Transmitter</b>	Yes
<b>Range</b>	-50°C to +150°C or 0°C to +400°C (either can be re-ranged if desired)
<b>Output</b>	4-20 mA output