

**Carcasa remota montada en la pared de 4-20 mA,  
entrada Pt100 con cable SSOB y sensor PRT -  
Longitud de cable de 2 m**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Carcasa moldeada en ABS de uso general, los puntos de fijacin externos proporcionan un medio sencillo y rpido de montar sensores directamente en los conductos, prescindiendo de la necesidad de placas de montaje separadas.

### **Necesita calibrar su sensor de temperatura?**

Trazable a los estndares de calibracin ISO17025 (UKAS), Labfacility ofrece la opcin de una calibracin trazable de 3 o 5 puntos utilizando nuestra instalacin de calibracin interna, puede seleccionar cualquier punto de temperatura entre -10 C y +200 C.

**Tenga en cuenta que si agrega una opcin de calibracin a su pedido, nos pondremos en contacto con usted con respecto a la informacin de entrega.**

### **Caractersticas tcnicas**

## **Specifications**

<b>Product Code</b>	XE-5663-001
<b>General Description</b>	Compatible with Pt100 platinum resistance thermometers
<b>Sensor Type</b>	Surface Mounting for remote installation
<b>Probe Diameter</b>	6mm
<b>Probe Length</b>	300mm
<b>Cable Length</b>	2m
<b>Cable Type</b>	Fiberglass stainless steel overbraid available
<b>Core / Strands</b>	7/0.2mm
<b>Head Type</b>	High Impact ABS Plastic Construction
<b>Thread</b>	Entry Threads M16
<b>Transmitter</b>	Yes
<b>Range</b>	0°C to +400°C
<b>Output</b>	4-20 mA output