

## Termmetro digital de catering tipo T



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

### Termmetro de catering

Termmetro de chef de alta precisin

Termopar tipo T

Cumple con HACCP,

-100 hasta + 280C

Funcin de retencin,

Termmetro impermeable IP67

Pantalla LCD dual

Fcil de usar: operacin con dos botones

Pantalla dual

Capacidad de apagado automtico.

Funciones:

Autocalibracin

Indicacin de batera baja

Se suministra con 1 batera PP3 de 9V

Rango de medicin -100 a +280 C

Alta precisin +/- 0.1% de lectura

Resolucin de +/-0,2 C, 0,1 C

El CA2005 ofrece una funcin de retencin y cuando se presiona este botn, la temperatura actual se mantendr y se transferir a la parte superior derecha de la pantalla. La pantalla principal continuar mostrando la temperatura en tiempo real.

### Caractersticas tcnicas

## Specifications

<b>Product Code</b>	XE-4090-001
<b>General Description</b>	High accuracy chef thermometer
<b>Type</b>	IP67 waterproof thermometer

<b>Accuracy</b>	+/-0.1% of Reading
<b>Resolution</b>	+/-0.2°C, 0.1°C Resolution
<b>Input / Ranges</b>	-100 to + 280°C
<b>Power Supply</b>	Supplied with 1 9V PP3 battery
<b>Sensor Type</b>	T