

## Termómetro de Legionella con temporizador de temperatura de agua integral



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

MM2008 Termómetro de entrada única con temporizador integral de 1 minuto y 2 minutos.

Este termómetro se desarrolla principalmente para su uso en la gestión de riesgos de legionela y ofrece la seguridad de que se toma la temperatura correcta al monitorear las temperaturas del agua caliente y fría.

### **Funciones**

Termopar de entrada única

Temporizador Integral para cuenta regresiva de 1 y 2 minutos

Rango de medición: -200 a +1,372 oC (dependiente del tipo de sonda)

Precisión: - - 0,1% de lectura a 0,2 oC

C/F Seleccionable

Revestimiento impermeable IP67

Indicación de batería baja

Tipos de termopar seleccionables: K & T

Autocalibración

Bota de goma incluida para una protección adicional contra danos

Incluye 1 x batería PP3

Requiere: Una sonda/sensor de termopar

### **Especificaciones**

## **Specifications**

<b>Product Code</b>	XE-4091-001
<b>General Description</b>	Measurement Range: -200 to +1,372°C (probe type dependant)
<b>Type</b>	Single input Thermometer with an integral 1 minute & 2 minute timer
<b>Accuracy</b>	± 0.1% of Reading ± 0.2°C
<b>Readings</b>	C/F Selectable
<b>Power Consumption</b>	Low Battery Indication
<b>Power Supply</b>	Includes 1 x PP3 Battery



# LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

<b>Standards Met</b>	IP67 Waterproof Casing
<b>Thermocouple Type</b>	K & T
<b>Thermometer</b>	Self-Calibrating
<b>Termination Type</b>	Requires a Thermocouple Probe/Sensor
<b>Accessories</b>	Rubber Boot included for extra protection against damage