


## RTD Plomo Rizado Rizado



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

## Specifications

<b>General Description</b>	RTD Retractable Curly Lead
<b>Cable Length</b>	2 or 5m (when fully extended)
<b>Cable Type</b>	RTD
<b>Core / Strands</b>	2 or 4 strand available
<b>Colour</b>	Black
<b>Insulation Material</b>	Polyurethane outer
<b>Termination Type</b>	Terminated with bare tails