


RTD Plomo Rizado Rizado



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Specifications

General Description	RTD Retractable Curly Lead
Cable Length	2 or 5m (when fully extended)
Cable Type	RTD
Core / Strands	2 or 4 strand available
Colour	Black
Insulation Material	Polyurethane outer
Termination Type	Terminated with bare tails