


## Enchufe conector de termopar miniatura FMTC-CU-M Tipo Cu



**Miniature in-line Plug**  
(dimensions in mm)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Enchufe conector de termopar miniatura FMTC-CU-M Tipo CU**

Con pines planos para adaptarse al tipo de termopar e ideales para uso general, todos los contactos estn polarizados para garantizar una conexin correcta. Estos conectores aceptarn cables de termopar de hasta 4 mm de dimetro.

**Caractersticas tcnicas****Specifications**

<b>Product Code</b>	XE-1370-001
<b>General Description</b>	Ideal for general purpose use. Contacts are polarised to ensure correct connection.
<b>Thermocouple Type</b>	IEC, ANSI, JIS Type CU
<b>Connector Type</b>	Plug
<b>Connector Size</b>	Miniature
<b>Colour</b>	White
<b>Dimensions</b>	See Schematic Drawing (Product Literature Tab)
<b>Max. Temperature</b>	220°C
<b>Pins</b>	Flat Pins
<b>+Positive Contact</b>	Copper
<b>-Negative Contact</b>	Copper
<b>Standards Met</b>	BS EN 50212:1997. RoHS