

## Sensor de temperatura RTD con conector Hirschmann - Longitud de la sonda de 75 mm



**'HIRSCHMANN PT100 PROBE ASSEMBLY'**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Un conjunto de sensor Pt100 rentable de 6 mm de diámetro alojado en un cuerpo de acero inoxidable 316 con conexión a proceso BSP de 1/2" y conector DIN instalado con cable de conexión.

**Características técnicas****Specifications**

<b>Product Code</b>	XE-5682-001
<b>General Description</b>	RTD Temperature Sensor & Hirschmann Connector
<b>Sensor Type</b>	PT100 Class A
<b>Minimum Span</b>	25°C *1
<b>Process Connection</b>	1/2" BSP
<b>Probe Diameter</b>	6.0mm
<b>Probe Length</b>	75mm
<b>Cable Length</b>	1m
<b>Cable Type</b>	PVC Insulated
<b>Max. Temperature</b>	250°C
<b>Min. Temperature</b>	-50°C
<b>IP Rating</b>	IP65
<b>Connection Type</b>	DIN connector Screw terminal
<b>Number of Wires</b>	4 wire connection