


Sensor de temperatura RTD con conector Hirschmann - Longitud de la sonda de 75 mm



'HIRSCHMANN PT100 PROBE ASSEMBLY'



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Un conjunto de sensor Pt100 rentable de 6 mm de diámetro alojado en un cuerpo de acero inoxidable 316 con conexión a proceso BSP de 1/2" y conector DIN instalado con cable de conexión.

Características técnicas**Specifications**

Product Code	XE-5682-001
General Description	RTD Temperature Sensor & Hirschmann Connector
Sensor Type	PT100 Class A
Minimum Span	25°C *1
Process Connection	1/2" BSP
Probe Diameter	6.0mm
Probe Length	75mm
Cable Length	1m
Cable Type	PVC Insulated
Max. Temperature	250°C
Min. Temperature	-50°C
IP Rating	IP65
Connection Type	DIN connector Screw terminal
Number of Wires	4 wire connection