


Conector de termopar en miniatura Conector en IAfÅ-nea IM-E-FS Tipo E IEC



Miniature in-line Socket
(dimensions in mm)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Conector de termopar en miniatura Conector en línea IM-E-FS Tipo E IEC

Ideales para uso general, todos los contactos están polarizados para garantizar una correcta conexión. Estos conectores aceptarán un cable termopar de hasta 4 mm de diámetro.

Especificaciones**Specifications**

Product Code	XE-1332-001
General Description	Ideal for general purpose use. Contacts are polarised to ensure correct connection.
Thermocouple Type	E IEC
Connector Type	Socket
Connector Size	Miniature
Colour	Violet
Max. Temperature	220°C
+Positive Contact	Nickel Chromium
-Negative Contact	Copper Nickel
Standards Met	BS EN 50212:1997. RoHS