

## Conector de cable rpido en miniatura Termopar Socket IM-K-FQ Tipo K IEC



**Miniature in-line Socket - Quick Wire**  
(dimensions in mm)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Conector de cable rpido en miniatura Termopar Socket IM-K-FQ Tipo K IEC**

Fcil conexin de jab-in, simplemente empuje el alambre y apriete el tornillo. Los contactos estn polarizados para evitar una terminacin incorrecta.

**Caractersticas tcnicas****Specifications**

|                            |  |
|----------------------------|--|
| <b>Product Code</b>        | XE-1112-001  |
| <b>General Description</b> | Easy jab-in connection, just push in wire & tighten screw. |
| <b>Thermocouple Type</b>   | K IEC  |
| <b>Connector Type</b>      | Quick Wire Socket  |
| <b>Connector Size</b>      | Miniature  |
| <b>Colour</b>              | Green  |
| <b>Max. Temperature</b>    | 220°C  |
| <b>+Positive Contact</b>   | Nickel Chromium  |
| <b>-Negative Contact</b>   | Nickel Alloy   |
| <b>Standards Met</b>       | BS EN 50212:1997. RoHS                                     |