


Bayoneta ajustable, cable trenzado de acero inoxidable de fibra de vidrio - tipo PRT - 2 m de longitud



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Bayoneta ajustable con resorte de 200 mm con tapa de fijacin ajustable (11,5 mm de diametro interno ranura nica).

Caracterstias tcnicas**Specifications**

| | |
|----------------------------|--|
| Product Code | XE-5655-001 |
| General Description | PRT with adjustable insertion depth |
| Sensor Type | Pt100 4-wire |
| Cable Length | 2M |
| Cable Type | Glassfibre insulated lead with stainless steel Overbraid |
| Core / Strands | 4 core 7/0.2mm |
| Cable Termination | Bare tails |
| Max. Temperature | 350°C |
| Min. Temperature | -50°C |
| Standards Met | IEC-60751 Class B |
| Bayonet cap | 11.5mm I.D. single slot |
| Spring | 6.0mm diameter x 200mm long |