


## Bayoneta ajustable Pt100, cable sobrebrazado de acero inoxidable Glassfibre - Tipo RTD / PRT



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Bayoneta ajustable PRT.**

Bayoneta ajustable (Pt100) PRT / RTD con muelle de 200 mm con tapa de fijación ajustable (11,5 mm I.D. ranura única)

Disponible en longitudes 1M o 2M.

**Specifications**

<b>General Description</b>	PRT with adjustable insertion depth
<b>Sensor Type</b>	Pt100 4 wire
<b>Cable Length</b>	1 or 2 metres
<b>Cable Type</b>	Glassfibre insulated lead with stainless steel Overbraid
<b>Core / Strands</b>	4 core 7/0.2mm
<b>Termination Type</b>	Bare Tails
<b>Max. Temperature</b>	+350°C
<b>Min. Temperature</b>	-50°C
<b>Standards Met</b>	IEC-60751 Class B
<b>Bayonet cap</b>	11.5mm I.D. single slot
<b>Spring</b>	6.0mm diameter x 200mm long