


Bayoneta ajustable Pt100, cable sobrebrazado de acero inoxidable Glassfibre - Tipo RTD / PRT



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Bayoneta ajustable PRT.

Bayoneta ajustable (Pt100) PRT / RTD con muelle de 200 mm con tapa de fijación ajustable (11,5 mm I.D. ranura única)

Disponible en longitudes 1M o 2M.

Specifications

General Description	PRT with adjustable insertion depth
Sensor Type	Pt100 4 wire
Cable Length	1 or 2 metres
Cable Type	Glassfibre insulated lead with stainless steel Overbraid
Core / Strands	4 core 7/0.2mm
Termination Type	Bare Tails
Max. Temperature	+350°C
Min. Temperature	-50°C
Standards Met	IEC-60751 Class B
Bayonet cap	11.5mm I.D. single slot
Spring	6.0mm diameter x 200mm long