


Conector de cable rpido en miniatura Zcalo de termopar FMTC-CU-FQ Tipo Cu



Miniature in-line Socket - Quick Wire
(dimensions in mm)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Conector de cable rpido en miniatura Termopar Socket FMTC-CU-FQ Tipo CU

Fcil conexin de jab-in, simplemente empuje el alambre y apriete el tornillo. Los contactos estn polarizados para evitar una terminacin incorrecta.

Caractersticas tcnicas**Specifications**

Product Code	XE-1382-001
General Description	Easy jab-in connection, just push in wire & tighten screw.
Thermocouple Type	IEC, ANSI, JIS Type CU
Connector Type	Quick Wire Socket
Connector Size	Miniature
Colour	White
Max. Temperature	220°C
+Positive Contact	Copper
-Negative Contact	Copper
Standards Met	BS EN 50212:1997. RoHS