




Sensores de temperatura de pastillas de freno de termopar IEC tipo K, J



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Estos sensores de termopar se utilizan para controlar la temperatura de las pastillas de freno hasta +400 C. Estn diseñados para ser incrustados en las almohadillas.

Cada sensor tiene 1 m de cable aislado de fibra de vidrio, con trenza de acero inoxidable. Disponible en termopar tipo K o J y con una seleccin de dimetros de punta. Cada punta est sellada con una pasta cermica de alta temperatura (Fortafix) que puede soportar temperaturas de hasta +400 C.

Specifications

General Description	Designed for monitoring the temperature of brake pads
Sensor Type	Grounded Junction
Thermocouple Type	K or J
Probe Diameter	3.2mm, 4.0mm or 6.0mm
Probe Length	Tip Length: 4mm
Cable Length	1m
Cable Type	Fibre glass with stainless steel Overbraid 7/0.2mm
Termination Type	Bare tails
Max. Temperature	+400°C
Features #1	Potting sealant: Ceramic paste (Fortafix)