

**Pt100 Termometro de resistencia dplex de 2 x 3  
hilos clase A, cabezal KNE - Longitud de la sonda  
de 150 mm**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Características técnicas****Specifications**

<b>Product Code</b>	XE-3656-001
<b>Sensor Type</b>	2 x Pt100 (100 Ohms @ 0°C) to IEC 751, Class A
<b>Probe Diameter</b>	6mm
<b>Probe Length</b>	150mm
<b>Sheath Material</b>	316 stainless steel
<b>Head Type</b>	KNE
<b>Termination Type</b>	6 way connection block for KNE head
<b>Max. Temperature</b>	+250°C (connection head @ 170°C)
<b>Min. Temperature</b>	-50°C
<b>IP Rating</b>	IP67
<b>Accuracy</b>	Class A
<b>Number of Wires</b>	2 x 3 Wire