



# LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

## Pt100 Duplex 2 x 3 clase de alambre A Termómetro de Resistencia, Cabeza KNE



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

## Specifications

<b>Sensor Type</b>	2 x Pt100 (100 Ohms @ 0°C) to IEC 751, Class A
<b>Probe Diameter</b>	6mm
<b>Probe Length</b>	100mm, 150mm, 200mm or 250mm
<b>Sheath Material</b>	316 stainless steel
<b>Head Type</b>	KNE
<b>Termination Type</b>	6 way connection block for KNE head
<b>Max. Temperature</b>	+250°C (connection head @ 170°C)
<b>Min. Temperature</b>	-50°C
<b>IP Rating</b>	IP67
<b>Accuracy</b>	Class A
<b>Number of Wires</b>	2 x 3 Wire