

## Termopar con aislamiento mineral de estilo industrial con cabezal de terminal KNE - Inconel tipo K (6 mm x 400 mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

## Specifications

|                              |  |
|------------------------------|--|
| <b>Product Code</b>          | XE-3248-001  |
| <b>General Description</b>   | Highly flexible, sheath can be bent/formed to suit many applications and processes |
| <b>Sensor Type</b>           | Mineral insulated Thermocouple fitted with KNE terminal head                       |
| <b>Thermocouple Type</b>     | Type K   |
| <b>Thermocouple Junction</b> | Insulated hot junction   |
| <b>Probe Diameter</b>        | 6.0mm diameter   |
| <b>Probe Length</b>          | 400mm Length   |
| <b>Sheath Material</b>       | Inconel  |
| <b>Head Type</b>             | KNE, Cast aluminium finish, screw type lid with chain                              |
| <b>Gland</b>                 | M16 x 1.5mm Plated brass gland fitted  |
| <b>Max. Temperature</b>      | +1100°C  |
| <b>Min. Temperature</b>      | -40°C  |
| <b>IP Rating</b>             | IP67   |