


Conector de fascia de termopar rectangular en miniatura Zcalo de fascia FMTC-CU-FF Tipo Cu



Miniature Fascia Socket (with nylon mounting clip)
(dimensions in mm)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Conector de termopar rectangular en miniatura Zcalo de fascia FMTC-CU-FF Tipo CU

Clip de nylon para una fijacin rpida y sencilla al panel. Todos los zcalos de montaje de la fascia estn polarizados por medio de pasadores de tamao desigual

Caractersticas tcnicas**Specifications**

Product Code	XE-1377-001
General Description	All fascia mounting sockets are polarised by means of unequally sized pins
Thermocouple Type	IEC, ANSI, JIS Type CU
Connector Type	Fascia mounting socket
Connector Size	Miniature
Colour	White
Max. Temperature	120°C
+Positive Contact	Copper
-Negative Contact	Copper
Standards Met	BS EN 50212:1997. RoHS