





## Conector de termopar en miniatura AM-K-M + F tipo K ANSI



**Miniature in-line Plug**  
(dimensions in mm)



**Miniature in-line Socket**  
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Conector de termopar en miniatura y enchufe AM-K-M + F tipo K ANSI**

Con clavijas planas para adaptarse al tipo de termopar e ideales para uso general, todos los contactos estn polarizados para garantizar una conexin correcta. Estos conectores aceptarn cables de termopar de hasta 4 mm de dimetro.

**Especificaciones****Specifications**

<b>Product Code</b>	XE-1135-001
<b>General Description</b>	Ideal for general purpose use. Contacts are polarised to ensure correct connection
<b>Thermocouple Type</b>	K ANSI
<b>Connector Size</b>	Miniature
<b>Connector Type</b>	Plug & Socket
<b>Colour</b>	Yellow
<b>Max. Temperature</b>	220°C
<b>Pins</b>	Flat Pins
<b>+Positive Contact</b>	Nickel Chromium
<b>-Negative Contact</b>	Nickel Alloy
<b>Standards Met</b>	BS EN 50212:1997. RoHS