



Conector de termopar estándar Enchufe y enchufe AS-B-M+F Tipo B ANSI



Standard in-line Plug
(dimensions in mm)



Standard in-line Socket
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Conector de termopar estándar Enchufe y enchufe AS-B-M+F Tipo B ANSI

Los contactos están hechos de verdaderas aleaciones de termopar. Adecuado para aplicaciones industriales

Especificaciones**Specifications**

Product Code	XE-1063-001
General Description	Contacts are polarised to prevent incorrect connection
Thermocouple Type	B ANSI
Connector Size	Standard
Connector Type	Plug & Socket
Colour	Grey
Max. Temperature	220°C
Pins	Solid Round Pins
+Positive Contact	Copper
-Negative Contact	Copper
Standards Met	BS EN 50212:1997. RoHS