




Enchufe estndar del conector del termopar IS-K-M Tipo K IEC



Standard in-line Plug
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Enchufe estndar del conector del termopar IS-K-M Tipo K IEC

Los contactos estn hechos de verdaderas aleaciones de termopar. Adecuado para aplicaciones industriales

Caractersticas tcnicas**Specifications**

| | |
|----------------------------|--|
| Product Code | XE-1400-001 |
| General Description | Contacts are polarised to prevent incorrect connection |
| Thermocouple Type | K IEC |
| Connector Type | Plug |
| Connector Size | Standard |
| Colour | Green |
| Max. Temperature | 220°C |
| Pins | Solid Round Pins |
| +Positive Contact | Nickel Chromium |
| -Negative Contact | Nickel Alloy |
| Standards Met | BS EN 50212:1997. RoHS |