


Conector de termopar estÃ¡fÃ¡ndar enchufe IS-E-M Tipo E IEC



Standard in-line Plug
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Conector de termopar estándar enchufe IS-E-M Tipo E IEC

Los contactos están hechos de verdaderas aleaciones de termopar. Adecuado para aplicaciones industriales

Especificaciones**Specifications**

Product Code	XE-1630-001
General Description	Contacts are polarised to prevent incorrect connection
Thermocouple Type	E IEC
Connector Type	Plug
Connector Size	Standard
Colour	Violet
Max. Temperature	220°C
Pins	Solid Round Pins
+Positive Contact	Nickel Chromium
-Negative Contact	Copper Nickel
Standards Met	BS EN 50212:1997. RoHS