



## Parche de caucho de silicona Termopar PFA Twin Twisted con terminacin de enchufe en miniatura - Tipo K JIS - 1m



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Caractersticas tcnicas****Specifications**

<b>Product Code</b>	XE-0555-001
<b>General Description</b>	Resistant to many oils & chemicals
<b>Sensor Type</b>	Thermocouple junction encapsulated in a silicone rubber patch
<b>Thermocouple Type</b>	K JIS
<b>Cable Length</b>	1m
<b>Cable Type</b>	PFA insulated twin-twist wire 1/0.2mm
<b>Termination Type</b>	Miniature flat pin plug
<b>Max. Temperature</b>	+150°C
<b>Min. Temperature</b>	-50°C
<b>Standards Met</b>	JIS (Japanese Industrial Standards)
<b>Colour</b>	Positive: Red / Negative: White
<b>Dimensions</b>	Length 40mm x Width 13mm x Height 5mm Rubber patch