


Parche de caucho de silicona Termopar PFA Twin Twisted con terminacin de enchufe en miniatura - Tipo K JIS - 5m



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Caractersticas tcnicas**Specifications**

Product Code	XE-0557-001
General Description	Resistant to many oils & chemicals
Sensor Type	Thermocouple junction encapsulated in a silicone rubber patch
Thermocouple Type	K JIS
Cable Length	5m
Cable Type	PFA insulated twin-twist wire 1/0.2mm
Termination Type	Miniature flat pin plug
Max. Temperature	+150°C
Min. Temperature	-50°C
Standards Met	JIS (Japanese Industrial Standards)
Colour	Positive: Red / Negative: White
Dimensions	Length 40mm x Width 13mm x Height 5mm Rubber patch