






# LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

## Termopar de parche de caucho de silicona JIS, PFA doble retorcido con tapn en miniatura o colas desnudas - tipo K



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Caractersticas tcnicas****Specifications**

<b>General Description</b>	Resistant to many oils & chemicals
<b>Sensor Type</b>	Thermocouple junction encapsulated in a silicone rubber patch
<b>Thermocouple Type</b>	K
<b>Cable Length</b>	1m, 2m, 5m
<b>Cable Type</b>	PFA insulated twin-twist wire 1/0.2mm
<b>Termination Type</b>	Miniature flat pin plug or Bare tails
<b>Max. Temperature</b>	+150°C
<b>Min. Temperature</b>	-50°C
<b>Standards Met</b>	JIS (Japanese Industrial Standards)
<b>Colour</b>	Positive: Red / Negative: White
<b>Dimensions</b>	Length 40mm x Width 13mm x Height 5mm Rubber patch