



# LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

## Termopar de disco de alambre fino, PFA con terminacin de cola desnuda - Tipo K JIS - 5m



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Caractersticas tcnicas****Specifications**

|                              |  |
|------------------------------|--|
| <b>Product Code</b>          | XE-0566-001  |
| <b>General Description</b>   | Exposed junction thermocouple with Copper Disc brazed at tip   |
| <b>Sensor Type</b>           | Copper disc, 6.35mm diameter x 0.3mm thick   |
| <b>Thermocouple Type</b>     | K JIS  |
| <b>Thermocouple Junction</b> | Grounded measuring junction  |
| <b>Cable Length</b>          | 5m   |
| <b>Cable Type</b>            | 1/0.2mm PFA Teflon® insulated twisted wire   |
| <b>Max. Temperature</b>      | 200°C  |
| <b>Min. Temperature</b>      | 0°C  |
| <b>Standards Met</b>         | JIS (Japanese Industrial Standards)  |
| <b>Colour</b>                | Positive: Red / Negative: White  |
| <b>Accessories</b>           | Supplied with high temperature adhesive kapton strips (13x18mm x 5) for quick and easy attachment to surface/process |