


Termopar de parche de caucho de silicona ANSI, PFA doble retorcido con tapn en miniatura o colas desnudas - Tipo K



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Características técnicas**Specifications**

General Description	Resistant to many oils & chemicals
Sensor Type	Thermocouple junction encapsulated in a silicone rubber patch
Thermocouple Type	Type K
Cable Length	1m, 2m, 5m
Cable Type	PFA insulated twin-twist wire 1/0.2mm
Termination Type	Miniature flat pin plug or Bare tails
Max. Temperature	+150°C
Min. Temperature	-50°C
Standards Met	ANSI MC96.1
Colour	Positive: Red / Negative: Yellow
Dimensions	Length 40mm x Width 13mm x Height 5mm Rubber patch