

**Termopar de pernos, acero inoxidable de fibra de vidrio sobre cable trenzado - JIS tipo K- Paso de rosca M8x1 mm - Tamao de perno de 25 mm**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Los sensores de pernos generalmente se atornillan en la superficie, lo que normalmente aumenta la velocidad de respuesta y la sensibilidad de la medicion. Apto para su uso hasta +250C.

**Caractersticas tcnicas****Specifications**

<b>Product Code</b>	XE-0591-001
<b>Thermocouple Type</b>	K JIS
<b>Thermocouple Junction</b>	Bolts have grounded junctions
<b>Cable Length</b>	2m
<b>Cable Type</b>	Glassfibre stainless steel over braided cable
<b>Cable Termination</b>	Bare Tails
<b>Max. Temperature</b>	250°C
<b>Standards Met</b>	JIS (Japanese Industrial Standards)
<b>Bolt thread/pitch</b>	M8 x 1mm
<b>Colour</b>	Positive: Red / Negative: White
<b>Length</b>	25mm