

Termopar de pernos, acero inoxidable de fibra de vidrio sobre cable trenzado - JIS tipo K- Paso de rosca M12x1 mm - Tamao de perno de 25 mm



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Los sensores de pernos generalmente se atornillan en la superficie, lo que normalmente aumenta la velocidad de respuesta y la sensibilidad de la medicion. Apto para su uso hasta +250C.

Caractersticas tcnicas**Specifications**

Product Code	XE-0595-001
Thermocouple Type	K JIS
Thermocouple Junction	Bolts have grounded junctions
Cable Length	2m
Cable Type	Glassfibre stainless steel over braided cable
Cable Termination	Bare Tails
Max. Temperature	250°C
Standards Met	JIS (Japanese Industrial Standards)
Bolt thread/pitch	M12 x 1mm
Colour	Positive: Red / Negative: White
Length	25mm