

**Termopar magnético de botón JIS tipo K de 1,9 kg,  
cable aislado con PFA con trenza de acero  
inoxidable que termina en enchufe en miniatura -**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Caractersticas tcnicas****Specifications**

<b>Product Code</b>	XE-0624-001
<b>General Description</b>	1.9kg Pull Button Magnet Thermocouple
<b>Sensor Type</b>	K JIS
<b>Magnet</b>	19.1mm Diameter / 12.7mm height - Max pull: 1.9 KG
<b>Cable Length</b>	1m
<b>Cable Type</b>	PFA Insulated Cable with stainless steel overbraid 7/0.2mm
<b>Termination Type</b>	Miniature Plug
<b>Max. Temperature</b>	250°C
<b>Min. Temperature</b>	-50°C
<b>Standards Met</b>	JIS (Japanese Industrial Standards)