

**Termopar magnético de botón de extracción tipo K  
JIS de 1,9 kg, cable aislado de PFA con trenza de  
acero inoxidable que termina en un enchufe en  
miniatura - 2 m**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Caracterstias tcnicas****Specifications**

<b>Product Code</b>	XE-0625-001
<b>General Description</b>	1.9kg Pull Button Magnet Thermocouple
<b>Thermocouple Type</b>	K JIS
<b>Magnet</b>	19.1mm Diameter / 12.7mm height - Max pull: 1.9 KG
<b>Cable Length</b>	2m
<b>Cable Type</b>	PFA Insulated Cable with stainless steel overbraid 7/0.2mm
<b>Termination Type</b>	Miniature Plug
<b>Max. Temperature</b>	250°C
<b>Min. Temperature</b>	-50°C
<b>Standards Met</b>	JIS (Japanese Industrial Standards)