



# LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

## Termopar de clip de cocodrilo ANSI con cable trenzado de acero inoxidable de fibra de vidrio - Tipo K



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Caractersticas tcnicas****Specifications**

<b>General Description</b>	Crocodile Clip Thermocouple with Glassfibre Stainless Steel Overbraided Cable
<b>Thermocouple Type</b>	K
<b>Thermocouple Junction</b>	Grounded measuring junction
<b>Cable Length</b>	1 or 3 metre
<b>Cable Type</b>	stainless steel over-braided PFA Teflon®
<b>Cable Termination</b>	Miniature Plug
<b>Max. Temperature</b>	250°C
<b>Min. Temperature</b>	-50°C
<b>Standards Met</b>	ANSI MC96.1