

## Termopar Crocodile Clip con cable trenzado de acero inoxidable de fibra de vidrio - JIS Tipo K - 3 m de longitud



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Caractersticas tcnicas****Specifications**

<b>Product Code</b>	XE-0631-001
<b>General Description</b>	Crocodile Clip Thermocouple with Glassfibre Stainless Steel Overbraided Cable
<b>Thermocouple Type</b>	K JIS
<b>Cable Length</b>	3m
<b>Cable Type</b>	Stainless steel over-braided PFA Teflon®
<b>Cable Termination</b>	Miniature Plug
<b>Max. Temperature</b>	250°C
<b>Min. Temperature</b>	-50°C
<b>Standards Met</b>	JIS (Japanese Industrial Standards)