

Termopar de clip de cocodrilo JIS con cable trenzado de acero inoxidable de fibra de vidrio - Tipo K



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Caractersticas tcnicas**Specifications**

General Description	Crocodile Clip Thermocouple with Glassfibre Stainless Steel Overbraided Cable
Thermocouple Type	K
Thermocouple Junction	Grounded measuring junction
Cable Length	1 or 3 metre
Cable Type	Stainless steel over-braided PFA Teflon®
Cable Termination	Miniature Plug
Max. Temperature	+250°C
Min. Temperature	-50°C
Standards Met	JIS (Japanese Industrial Standards)