

**Termopar conectado a tierra fabricado en tubo de  
acero inoxidable, cable trenzado de acero  
inoxidable de fibra de vidrio - Tipo K JIS - 4,76  
mm de diámetro y 25 mm de longitud**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Caractersticas tcnicas****Specifications**

<b>Product Code</b>	XE-0650-001
<b>Sensor Type</b>	316 stainless steel tube in various sizes
<b>Thermocouple Type</b>	K JIS
<b>Thermocouple Junction</b>	Grounded at tip for fast response
<b>Cable Length</b>	2m
<b>Cable Type</b>	Glassfibre stainless steel overbraided cable
<b>Cable Termination</b>	Bare Wire Tails
<b>Max. Temperature</b>	350°C
<b>Min. Temperature</b>	-60°C
<b>Standards Met</b>	JIS (Japanese Industrial Standards)
<b>Dimensions</b>	4.76mm Diameter x 25mm Length