

**ANSI Tipo K 4.5kg Pull Power (Herradura)
Termopar magnético, cable aislado PFA con trenza
de acero inoxidable que termina en miniatura o
enchufe estándar**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Características técnicas**Specifications**

General Description	4.5kg Pull Power (Horseshoe) Magnet Thermocouple
Sensor Type	Magnet Thermocouple
Thermocouple Type	K
Magnet	30 x 20 x 20mm (L x H x W). Max Pull: 4.5kg
Cable Length	1m or 2m
Cable Type	PFA Teflon® insulated with stainless steel over-braid
Core / Strands	7/0.2mm
Termination Type	Miniature or Standard Plug
Max. Temperature	+250°C
Min. Temperature	-50°C
Standards Met	ANSI MC96.1