

ANSI Tipo K 9kg Pull Power (herradura) Termopar magnético, cable aislado PFA con trenza de acero inoxidable que termina en miniatura o enchufe estándar



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Caractersticas tcnicas**Specifications**

General Description	Rapid, reliable temperature sensing on ferrous surfaces
Thermocouple Type	K
Thermocouple Junction	Grounded at tip
Magnet	25 x 25 x 40mm (HxWxL) with powerful 9kg pull
Cable Termination	1.5 or 3 metres PFA Teflon® insulated with stainless steel over-braid, fitted plug termination
Termination Type	Fitted Mini or Standard Plug
Max. Temperature	+250°C
Min. Temperature	-50°C
Standards Met	ANSI MC96.1