

**ANSI Tipo K 11.8kg Pull Power (herradura)
Termopar magnético, cable aislado PFA con trenza
de acero inoxidable que termina en miniatura o
enchufe estándar**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Características técnicas**Specifications**

General Description	11.8kg Pull Power (Horseshoe) Magnet Thermocouple
Sensor Type	Magnet Thermocouple
Thermocouple Type	K
Magnet	45 x 30 x 30mm (L x H x W) - Max pull: 11.8 KG
Cable Length	1.5m or 3m
Cable Type	PFA Insulated Cable with Stainless Steel Over-Braid
Core / Strands	7/0.2mm
Cable Termination	Miniature or Standard Plug
Max. Temperature	+250°C
Min. Temperature	-50°C
Standards Met	ANSI MC96.1